



Product Ordering Information

ON-Q® Custom Procedure Kits

Create your own Custom Procedure Kit with the ON-Q products of your choice—the ultimate choice for convenience and flexibility. Simply select the ON-Q pump, catheter(s), needle(s) and tunneler to design a custom kit that best meets your needs.

With Custom Procedure Kits, you also have the option to label the kits according to doctor, facility and/or pre-selected procedure. This enables your staff to easily identify kits used by specific doctors or quickly find the equipment needed for a specific procedure in one convenient package. Call your I-Flow representative for more information.

ON-Q® C-bloc®

Simple and reliable system for effective continuous nerve blocks—to get patients back to normal faster.

ON-Q C-bloc slowly infuses local anesthetic near a nerve for regional anesthesia and post-operative pain relief. ON-Q C-bloc with ONDEMAND™ has an added feature that allows patients to give a 5 ml bolus in addition to the continuous flow rate. The Select-A-Flow® feature allows you to adjust the flow rate according to your patient's needs.

Delivery Time* (in days)	Fill Volume x Flow Rate	Model	Per Package
3	400 ml x 5 ml/hr	CB001	5
2	400 ml x 10 ml/hr	CB002	5
ONDEMAND			
< 3	400 ml x 5 ml/hr with 5 ml bolus (60 minute refill)	CB003	5
Select-A-Flow			
Varies	400 ml x 2-14 ml/hr	CB004	5
Varies	600 ml x 2-14 ml/hr	CB6004	5
Varies	600 ml (dual; 1-7 ml/site)	CB6007	5
ONDEMAND & Select-A-Flow			
Varies	400 ml x 1-7 ml/hr with 5 ml bolus (60 minute refill)	CB005	5
Varies	400 ml x 2-14 ml/hr with 5 ml bolus (30 minute refill)	CB006	5

* Delivery time information is approximate and based on nominal specifications. Refer to Directions for Use for complete product information.


ON-Q® SilverDressing™

The thin film dressing with a powerful “one-two” punch of physical and antimicrobial barriers against infection.

First, the dressing creates a physical barrier that helps shield the wound site from infection. Second, the ionic silver is 99.99% effective in killing or inhibiting the growth of a wide range of microorganisms.**

Description	Model	Per Package
Antimicrobial Dressing 4 in x 4.5 in	SD4X4.5	25
Antimicrobial Dressing 2 in x 3 in	SD2X3	25

** ON-Q SilverDressing is not intended as a treatment for clinical infections.



ON-Q® PainBuster®

Significantly better, non-narcotic relief that's Redefining Recovery™.

ON-Q PAINBUSTER KITS



Delivering local anesthetic automatically and continuously to the surgical site for up to 5 days—providing effective pain relief.

Delivery Time* (in days)	Fill Volume x Flow Rate	Standard Catheter	Soaker® Catheter 1 in (2.5 cm)	Soaker Catheter 2.5 in (6.5 cm)	Soaker Catheter 5 in (12.5 cm)	Soaker Catheter 10 in (25 cm)	Per Package
5	65 ml x 0.5 ml/hr	PM001	PM031	PM011	—	—	5
2	100 ml x 2 ml/hr	PM002	PM032	PM012	PM022	—	5
5	270 ml x 2 ml/hr	PM003	—	PM013	PM023	—	5
2¼	270 ml x 5 ml/hr	PM004	—	PM014	PM024	—	5
3	270 ml x 4 ml/hr (dual; 2 ml/site)	PM005	PM035	PM015	PM025	PM045	5
4	100 ml x 1 ml/hr	PM006	PM036	—	—	—	5
3½	400 ml x 5 ml/hr	—	—	—	PM026	PM046	5
4	400 ml x 4 ml/hr (dual; 2 ml/site)	—	—	PM018	PM028	PM048	5
2	100 ml x 2 ml/hr (dual; 1 ml/site)	—	PM039	—	—	—	5

ON-Q PAINBUSTER KITS WITH ON-Q SILVERSOAKER™

The same ON-Q you know and trust with an added layer of protection for even better outcomes.

Delivery Time* (in days)	Fill Volume x Flow Rate	Soaker Catheter 1 in (2.5 cm)	Soaker Catheter 2.5 in (6.5 cm)	Soaker Catheter 5 in (12.5 cm)	Soaker Catheter 10 in (25 cm)	Per Package
2	100 ml x 2 ml/hr	PM032-A	PM012-A	PM022-A	—	5
5	270 ml x 2 ml/hr	—	PM013-A	PM023-A	—	5
2¼	270 ml x 5 ml/hr	—	PM014-A	PM024-A	—	5
3	270 ml x 4 ml/hr (dual; 2 ml/site)	—	PM015-A	PM025-A	PM045-A	5
4	100 ml x 1 ml/hr	PM036-A	—	—	—	5
3½	400 ml x 5 ml/hr	—	—	PM026-A	PM046-A	5
4	400 ml x 4 ml/hr (dual; 2 ml/site)	—	PM018-A	PM028-A	PM048-A	5



ON-Q PUMP-ONLY KITS (See ON-Q Expansion Kits below for corresponding Soaker® Catheters and needles)

Delivery Time* (in days)	Fill Volume x Flow Rate	Model	Per Package
5	65 ml x 0.5 ml/hr	P65X0.5	5
2	100 ml x 2 ml/hr	P100X2	5
5	270 ml x 2 ml/hr	P270X2	5
3	270 ml x 4 ml/hr	P270X4	5
2¼	270 ml x 5 ml/hr	P270X5	5
3	270 ml x 4 ml/hr (dual; 2 ml/site)	P270X4D	5
4	400 ml x 4 ml/hr	P400X4	5
3	400 ml x 5 ml/hr	P400X5	5
4	400 ml x 4 ml/hr (dual; 2 ml/site)	P400X4D	5
2	100 ml x 2 ml/hr (dual; 1 ml/site)	P100X2D	5
≤4	400 ml x 6 ml/hr (triple; 2 ml/site)	P400X6T	5
2–5	Select-A-Flow 270 ml x 1-7 ml/hr	SAF01	5

ON-Q ONDEMAND™ SYSTEM

Post-op pain relief delivered on demand.

Delivery Time* (in days)	Fill Volume x Flow Rate	Fenestrated Catheter (1 in)	Per Package
2–5	270 ml x 2 ml/hr (5 ml bolus) (60 minute refill)	PMB01	5

ON-Q SOAKER & ON-Q SILVERSOAKER CATHETER EXPANSION KITS

The patented Soaker Catheter provides an even distribution of local anesthetic for optimal post-op pain relief.

Description	Soaker Model	SilverSoaker Model	Per Package
Soaker 2.5 in (6.5 cm), T-Peel Needle 3.25 in	PM010	PM010-A	5
Soaker 5 in (12.5 cm), T-Peel Needle 6 in	PM020	PM020-A	5
Soaker 1in (2.5 cm), T-Peel Needle 3.25 in	PM030	PM030-A	5
Soaker 10 in (25 cm), T-Peel Needle 8 in	PM040	PM040-A	5
Soaker 7.5 in (19 cm), T-Peel Needle 8 in	—	PM050-A	5

ON-Q PAINBUSTER CK KITS

For surgical procedures with large incisions.

Description**	Soaker Model	SilverSoaker Model	Per Package
PM025, ACC03	—	CK025-A	5
PM025, PM040, ACC10S	CK045	CK045-A	5
PM028, PM040, ACC10S	CK048	CK048-A	5

** -A models include ON-Q SilverSoaker.





ON-Q TUNNELING SYSTEM

Taking effective pain relief a step beyond.

Description	Model	Per Package
Reusable Tuner 11 GA x 8 in	ACC07	1
Tuner Sheath 11 GA x 8 in	ACC07S	10
Reusable Tuner 17 GA x 8 in	ACC08	1
Tuner Sheath 17 GA x 8 in	ACC08S	10
Reusable Tuner 11 GA x 12 in	ACC09	1
Tuner Sheath 11 GA x 12 in	ACC09S	10
Reusable Tuner 16 GA x 12 in	ACC10	1
Tuner Sheath 16 GA x 12 in	ACC10S	10
Reusable Tuner 17 GA x 3.25 in	ACC12	1
Tuner Sheath 17 GA x 3.25 in	ACC12S	10
Reusable Tuner 17 GA x 5 in	ACC13	1
Tuner Sheath 17 GA x 5 in	ACC13S	10
Disposable Tuner & Sheaths 17 GA x 3.25 in	T17X3.25	5
Disposable Tuner & Sheaths 17 GA x 5 in	T17X5	5
Disposable Tuner & Sheaths 17 GA x 8 in	T17X8	5
Disposable Tuner & Sheaths 16 GA x 12 in	T16X12	5
Disposable Tuner & Sheaths 11 GA x 8 in	T11X8	5
Disposable Tuner & Sheaths 11 GA x 12 in	T11X12	5
Disposable Fluid Tuner & Sheaths 11 GA (12 in)	FT11X12	5
Disposable Fluid Tuner & Sheaths 11 GA (8 in)	FT11X8	5
Disposable Fluid Tuner & Sheaths 16 GA (12 in)	FT16X12	5
Disposable Fluid Tuner & Sheaths 16 GA (8 in)	FT16X8	5

ON-Q ACCESSORIES


Description	Model	Per Package
T-Peel Needle 3.25 in	ACC01	10
T-Peel Needle 6 in	ACC02	10
T-Peel Needle 8 in	ACC03	10
Large E-Clips (For use with 270 ml pumps.)	ACC04	25
Filling Extension Sets	ACC05	5
Reusable PowerFiller	ACC11	1
Locking Carrying Cases	ACC14	5



.....

Contact your ON-Q representative
for ordering information

.....



FOR MORE PRODUCT RELATED OR ORDERING INFORMATION CALL

800.448.3569 or **949.206.2700**

20202 Windrow Drive | Lake Forest | California 92630

*FOR A LIST OF SUPPORTING STUDIES VISIT

www.iflo.com

www.AskYourSurgeon.com