

FLUSH

Urokinase

63 Syringes  
for 7 Days

Now There is a Solution...



I-FLOW®

One • Step KVO

S.A.S.H.

occlusions

pole mounted pumps

VIALS

infections

SYRINGES



# One•Step KVO™

## Elastomeric pump

*Sometimes,*  
*in order to solve a very difficult problem,*  
*you have to think outside of the box.*  
*Whether you are looking for ways to*  
***reduce occlusions and infections,***  
***substantially reduce your flushing costs,***  
***or just free up nursing time...***  
*the One•Step KVO is the answer!*

### Potential Clinical Benefits


- Reduces incidence of catheter occlusions
- Decreases risk of infections
- Reduces risk of needle sticks

### Lower Costs

- Reduces nursing labor
- Reduces flushing costs up to 80%
- Eliminates associated costs of complications (Urokinase, catheter replacement, etc.)

### Patient Benefits

- Easily administered by patients or caregivers
- Completely ambulatory
- Ideal for pediatric patients



Description	Model No.	Units/Case	Days of Infusion	Flow Rate	Fill Volume
One•Step KVO	EY060005	24	4	0.5 ml/hr	60 ml
One•Step KVO	EY110005	24	7	0.5 ml/hr	110 ml
One•Step KVO w/Check Valve	EV060005	24	4	0.5 ml/hr	60 ml
One•Step KVO w/Check Valve	EV110005	24	7	0.5 ml/hr	110 ml
E-Clip Clothing Clip	5001087	100			
Carrying Case	1400808	1			

For Customer Service  
 Call: 1.800.448.3569  
 949.206.2700  
[www.i-flowcorp.com](http://www.i-flowcorp.com)  
 1302247A / 2/99

CE European Representative:  
 MPS Medical Product Service GmbH  
 0 1 2 3 Borggasse 20, 35619 Braunfels, Germany

A PRODUCT OF  
  
 I-FLOW CORPORATION  
 LAKE FOREST, CA 92630  
 U.S.A.