



Double Protection against Microbes

Introducing **ON-Q® SilverDressing™**. The thin film dressing with a powerful “one-two” punch of physical and antimicrobial barriers against microorganisms.*

First, the dressing creates a physical barrier that helps shield the wound site from microbes. Second, the ionic silver is **99.99% effective** in inhibiting the growth of a wide range of microorganisms including MRSA¹, for up to 7 days.

WHY SMART
CLINICIANS ARE CHOOSING
ON-Q SilverDressing

- Breathability means great comfort
- 99.99% antimicrobial effectiveness against a wide range of microorganisms¹
- A thin-film design you're used to, plus an added antimicrobial barrier

VISIT US ONLINE AT MyON-Q.com

ON-Q SilverDressing™
Redefining Recovery™

Silver is the new Gold

Knock out the high cost of fighting healthcare-associated infections (HAI).


375,000 reasons to invest in silver.

The CDC reports there are nearly 375,000 surgical site infections each year.² The cost of treating surgical site infections (SSI) can add thousands of dollars to hospital bills. Since October 1, 2008, Medicare no longer reimburses for preventable healthcare-associated infections (HAI).³ Now is the time to protect your patients and your facility!



ON-Q SilverDressing is strong yet flexible. It will effectively secure an ON-Q SilverSoaker catheter, and remain secure, even when the dressing is placed over a joint.

Double Protection Against Microbes



- A** Live microbes bouncing off the dressing may be killed
- B** Dead microbes on the skin around the wound
- C** Live microbes trying to get under the dressing may be killed

ON-Q SilverDressing's silver ions are released over time. It delivers extra protection against microbes—and the spread of microbes to others.

PERIPHERAL IV




PICC/MIDLINE




EPIDURAL



CENTRAL LINE



NOW available in two convenient sizes!



MODEL	SD2X3	SD4X4.5
PER CASE	25	25

ON-Q SilverDressing™
Redefining Recovery™

There are inherent risks in all medical devices, please refer to the product labeling for **Indications, Precautions, and Instructions for Use/Removal**. This medical device is not intended as a treatment for clinical infection. Product coloration may change but this does not affect the performance of the medical device.



FOR MORE PRODUCT RELATED INFORMATION CALL
800.448.3569 or 949.206.2700

20202 Windrow Drive | Lake Forest | California 92630

Rx only. I-Flow and ON-Q are registered trademarks of I-Flow Corporation. SilverDressing and Redefining Recovery are trademarks of I-Flow Corporation. ©2011 I-Flow Corporation. All rights reserved.

* ON-Q SilverDressing is not intended as a treatment for clinical infections.
1. Data maintained by I-Flow Corporation. Available on request. (n=3)
2. <http://www.cdc.gov/ncidod/dhqp/hai.html>
3. <http://www.nytimes.com/2007/08/19/washington/19hospital.html>